

CADPRO[®] Advantage With SmartFeeder[®]XL

X-Ray Film Digitizer for Mammography Priors, Digitized Film CAD

The VIDAR **CAD PRO Advantage** film digitizer was specifically designed and engineered for use in mammography CAD (computer-aided diagnosis) systems. CAD systems incorporate advanced pattern recognition and image analysis, requiring a film digitizer that meets the highest standards of image quality and reliability. By meeting these rigorous requirements, CAD PRO *Advantage* is an ideal fit for mammography facilities transitioning to digital technology. It provides the ability to quickly digitize prior studies and, when coupled with the appropriate software, uses the DICOM MG standard to identify the specific mammography view contained on each film. See VIDAR's Whitepaper "Realizing the Benefits of Full-Field Digital Mammography; *The Role of the Film Digitizer*" for more information.

In order to read a new mammography study, radiologists and mammographers may require the two most recent prior studies. The ongoing transition to digital technology in mammography often requires radiologists to view both film-based priors and new digital images for the same returning patient read. The CAD PRO *Advantage* is able to bridge the gap in the transition to digital mammography – allowing for the digitization of film-based prior studies so that radiologists can view these priors, with their specific hanging protocols, on the mammography reading workstation. Radiologists and mammographers have limited time for reading studies, and the need to move from digital workstation to light box and back is cumbersome and time-consuming – creating a less-than-ideal clinical and diagnostic environment.

Responding to Marketplace Needs With SmartFeeder[®] XL

With the introduction of the **SmartFeeder XL**, VIDAR proves its continued commitment to customer needs by extending the digitizer's capabilities beyond analog film. With the added capabilities of the SmartFeeder XL, both analog and printed mammography films can be digitized, ensuring a smooth workflow and increased productivity for the mammography team.

Technology Advantage and Innovation

The CAD PRO *Advantage* offers unmatched spatial and contrast resolution, allowing radiologists to be confident that the image data received from the film digitizer is an exact representation of the original film. VIDAR's Automatic Digitizer Calibration (ADC) prompts the film digitizer to calibrate automatically before every film digitized and user intervention is not needed to maintain image quality. The CAD PRO *Advantage* includes a closed-loop quality assurance system to ensure the digitizer sub-system meets all FDA/MQSA requirements. This combination ensures excellent grayscale reproduction.



CAD PRO Advantage with SmartFeeder XL

Meeting Customer Requirements

The CAD PRO *Advantage* includes a modular 50-sheet feeder. The advanced feeder design eliminates film jams, double feeds, and film pick-up problems that occur with other digitizers. Most importantly, the CAD PRO's feeder is designed for continuous case loading, allowing the user to add cases during the digitizing process.

The CAD PRO *Advantage* can digitize a four-film study in 93 seconds, 30% faster than previous versions, further increasing productivity. Reliability is ensured with the LED light source which improves repeatability of the device as well as longevity, and lowers the total cost of ownership. The modular design of the CAD PRO allows easy access for maintenance cleaning and light source field replacement, minimizing downtime.

CADPRO[®] Advantage

With SmartFeeder[®] XL

Film Size	Nominal Resolution	Pixels	Spot Size (µm)	DPI	Line pairs per mm	Digitizing Speed
24 cm x 18 cm (analog)	5K x 4K	5376 x 4032	44	570	11	17 seconds
24 cm x 30 cm (analog)	5K x 7K	5376 x 6720	44	570	11	29 seconds
10" x 8" (printed)	6K x 5K	5700 x 4560	44	570	11	20 seconds
10" x 12" (printed)	6K x 7K	5700 x 6870	44	570	11	29 seconds

Clinical Optical Density Range	.05 to 4.2
DMAX:	4.8
Bit Depth	32-bit mapped to 16-bit (65,536) or 12-bit (4,096) Grayscale Output
MTBF	≥50,000 hours
Film Sizes	18 cm x 24 cm, 24 cm x 30 cm, 8" x 10" and 10" x 12" Thickness: 0.006" to 0.008"
Auto Film Feeder	Continuous Loading SmartFeeder XL — Modular with 50-film capacity (mixed sized, no presorting necessary)
Translation Table	Linear OD
Geometric Accuracy	Better than 1% or 2 pixels, whichever is greater, in both axes
Scan Rate	233 lines/second (160 films per hour, 18cm x 24cm)
Hardware Interface	USB 2.0
Software	Windows [®] scanning modules available
Power Requirements	Voltage: 85~264 vac Frequency: 47~63 Hz Power: ≤100 watts
Operating Environment	50° to 95° F (10° to 35° C), 20% to 85% relative humidity non-condensing
Storage Environment	0° to 140° F (-18° to 60° C), 20% to 85% relative humidity non-condensing
Illuminator	LED Illuminator, > 500,000 scans
Detector	Solid-state, next-generation High Definition CCD (HD-CCD [®])
Dimensions	With Feeder & Exit Tray: 19" W x 21.25" D x 25.5" H (483 mm x 540 mm x 648 mm) Without Feeder & Exit Tray: 19" W x 14.25" D x 12.75" H (483 mm x 362 mm x 324 mm) Shipping: 24" W x 29" L x 24" H (610 mm x 737 mm x 610 mm)
Weight	47 lbs. (21 kg); shipping weight: 62 lbs (28 kg)

Specifications are subject to change without notice

365 Herndon Parkway
Herndon, VA USA 20170
www.filmdigitizer.com

Phone: +1.703.471.7070
Toll-free: 1.800.471.7226
Fax: +1.703.471.7665



VIDAR, CAD PRO, and SmartFeeder are registered trademarks of VIDAR Systems Corporation. All other product names are registered marks of their respective parent companies.