

Vidar

Digitizer Tools Version 2.3

Digitizer Diagnostics
Firmware Download
Passcode Editor

© Copyright 2013 by VIDAR Systems Corporation. All rights reserved.

No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written permission of VIDAR Systems Corporation.

VIDAR will provide support as per the terms and conditions set forth in the license agreement. Other terms and conditions may have been presented by VIDAR resellers. If these terms differ from those offered by VIDAR, it is the responsibility of the reseller, not VIDAR to fulfill these obligations. Further, support will only be given in the following situations:

- Installation and modifications are made by qualified personnel,
- The equipment is used in accordance with the instructions in this manual.
- No other software or hardware other than those sanctioned by VIDAR or distributed by VIDAR with Clinical Express are used on the Clinical Express workstation.

Contents

- How to use this manual 4**
- Digitizer Diagnostics..... 5**
 - Intended Use 5
 - Installing Digitizer Diagnostics 5
 - Using Digitizer Diagnostics 5
- Firmware Download 10**
 - Intended Use 10
 - Installing Firmware Download..... 10
 - Using Firmware Download..... 10
- Passcode Editor 13**
 - Intended Use 13
 - Installing Passcode Editor 13
 - Using Passcode Editor 13
- Troubleshooting 17**
- Appendix 24**

How to use this manual

This manual provides step-by-step instructions for installing and using these setup programs for Vidar film digitizers using **Windows XP, and Windows 7 (32 and 64 Bit)**:

- **Digitizer Diagnostics** comprehensively tests an attached digitizer.
- **Firmware Download** displays an attached digitizer's current firmware versions and assists in updating firmware.
- **Passcode Editor** enables a new passcode to be assigned to a digitizer.

NOTE: Note that all screen shots are used for instruction only. The versions may not reflect the current version of the tool being used, however the process is the same.

Digitizer Diagnostics

Intended Use

Digitizer Diagnostics performs comprehensive testing of an attached film digitizer. When testing is complete, you can send the results to Vidar Technical Support for analysis.


Installing Digitizer Diagnostics

1. Download *DigitizerDiagnostic.exe* from the Vidar Web site. www.filmdigitizer.com. Browse to the Support page and then browse to the Troubleshooting page. Save the file in a convenient location on the local computer.
2. Double-click *DigitizerDiagnostic.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the on screen instructions to complete the installation.

Using Digitizer Diagnostics

Note: The film digitizer must be powered on for at least 30 minutes before running diagnostic tests.

1. From the Desktop double click > **Digitizer Diagnostics**.
2. You will see the warning below. Click **OK** to continue.




3. Wait while the application loads.




Note: The program will notify you if it cannot communicate with the digitizer. In this case, go to the "Troubleshooting" chapter.

4. Click the **Start** button in the lower right corner to start the diagnostics.



5. Wait while the diagnostic tests are performed.



6. When the tests are complete, click the **Contact Form** button in the lower right corner.


7. Fill in the fields in the **Contact** form, then click **Next >>**.

The screenshot shows the Digitizer Diagnostics 2.2 application window. The title bar reads "Digitizer Diagnostics 2.2" with help and close buttons. The main window has a dark blue background. On the left, a vertical sidebar contains a progress indicator with three green checkmarks: "System check", "Requirements check", and "Run digitizer tasks", followed by the word "Done". The top right of the main area displays "Model: DiagnosticPRO Advantage (USA model)" and "Serial number: 361029". The central part of the window is a white form titled "Contact Information". Below the title, it says "Please fill out the form below and click the Next button below. Fields in bold are mandatory." The form contains the following fields: "Digitizer Serial No." (361029), "Contact Name" (Joe C Customer), "Company Name" (Customer), "Address" (Any Town USA), "Phone" (555-555-5555) with a note "(Include country and area codes)", "Email Address" (joe.customer@customer.com), and "Comments" (Unit will not Scan). At the bottom of the form are two buttons: "<< BACK" and "NEXT >>". In the bottom left corner of the application window is the Vidar Systems Corporation logo with the tagline "The Image of Reliability".


Digitizer Diagnostics generates a report of its findings. The report is located in the *My Documents\Film Digitizer\Report* folder. (A typical report is reproduced in the Appendix.)

8. Send the diagnostic report to Vidar using one of the following methods:
- If Digitizer Diagnostics detects e-mail capability on the computer, a new message window automatically appears with the diagnostic report attached. Send the message.
or
 - If Digitizer Diagnostics does not detect e-mail capability on the computer, you can copy the report file to a computer that has e-mail capability, then send the report as an attachment to medtech@vidar.com.
or
 - Print the report and fax it to 1.703.471.7665.

9. *To reset the digitizer:* Click **OK** in the dialog shown below. Digitizer Diagnostics will automatically close when the reset is complete or
To exit the application without resetting the digitizer: Click **Cancel**.



Note: You must reset the film digitizer before attempting to digitize films.



Firmware Download

Intended Use

Firmware Download determines a film digitizer's current firmware versions, and enables you to update firmware if newer versions are available.


NOTE: This should never run unless instructed to do so by a VIDAR authorized agent.

Installing Firmware Download

1. Download *FirmwareDownload.exe* from the Vidar FTP or PartnerNet sites. Save the file in a convenient location on the local computer.
2. Double-click *FirmwareDownload.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the on screen instructions to complete the installation.

Using Firmware Download

1. From the Desktop double click on **>Firmware Download**.
2. You will see the warning below. Click **OK** to continue.




3. Wait while the application loads.



Note: The program will notify you if it cannot communicate with the digitizer. In this case, go to the “Troubleshooting” chapter.


4. Firmware Download first reads and displays the digitizer’s firmware versions.

- If the **Status** column shows “valid” for all firmware versions, then no further action is required. Exit Firmware Download.
- If firmware updates recommended:
 - a. Open the drop-down list and select the desired firmware package to install.
 - b. Click the **Start** button in the lower right corner.



5. Wait while the firmware is updated.

CAUTION: Do not turn off the digitizer or computer during the firmware updating process.



6. When firmware updating is complete, the digitizer will be automatically reset. When resetting is complete, Firmware Download will automatically close.

Passcode Editor

Intended Use

Passcode Editor enables you to assign a new passcode to a film digitizer.


NOTE: This should never run unless instructed to do so by a VIDAR authorized agent.

Installing Passcode Editor

1. Download *PasscodeEditor.exe* from the Vidar FTP or PartnerNet sites. Save the file in a convenient location on the local computer.
2. Double-click *PasscodeEditor.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the onscreen instructions to complete the installation

Using Passcode Editor

1. **Start > All Programs > Vidar >Passcode Editor.**
2. You will see the warning below. Click **OK** to continue.



3. Wait while the application loads.



Note: The program will notify you if it cannot communicate with the digitizer. In this case, go to the “Troubleshooting” chapter.

4. In the **New password** field, enter the new password.

The screenshot shows the 'Passcode Editor 2.2' application window. The title bar includes a question mark icon and window control buttons. The main area displays the following information:

- Model: TeleradPRO Edge
- Serial number: 250004

On the left side, there is a vertical menu with the following items:

- System check (with a green checkmark)
- Requirements check (with a green checkmark)
- Run digitizer tasks
- Done

The main content area has a blue header with an information icon and the text 'Enter a valid passcode.'. Below this is a warning box with the following text:

Warning: Loading an incorrect passcode into the digitizer will damage the product with repair needed at VIDAR's headquarters. Only use this program if instructed to by VIDAR or your authorized reseller. This application is NOT capable of checking validity of the entire passcode. Double check the entry to avoid service issues. After the passcode is changed, the software will automatically reset the film digitizer. This is necessary and may take several minutes.

Below the warning box are three input fields:

- Old passcode: (A red arrow points to this field)
- New passcode:
- Confirm new passcode:

The VIDAR logo is visible in the bottom left corner of the application window.

The word “invalid” will appear if you enter an unacceptable passcode.

Old password	<input type="text" value="MKLIZAMBVLFLBN"/>	
New password	<input type="text" value="RAB901BNTW0291"/>	invalid ←
Confirm new password	<input type="text"/>	

5. Enter the new passcode again in the **Confirm new password** field. If the contents of the **Confirm new password** field match the contents of the **New password** field, the word “matched” will appear.

Passcode Editor 2.2

Model: TeleradPRO Edge
Serial number: 250004

System check ✓
Requirements check ✓
Run digitizer tasks
Done

Click the UPDATE button below to continue.

Warning: Loading an incorrect passcode into the digitizer will damage the product with repair needed at VIDAR's headquarters. Only use this program if instructed to by VIDAR or your authorized reseller. This application is NOT capable of checking validity of the entire passcode. Double check the entry to avoid service issues. After the passcode is changed, the software will automatically reset the film digitizer. This is necessary and may take several minutes.

Old passcode

New passcode valid


Confirm new passcode matched

UPDATE


The word “unmatched” will appear if the contents of the **Confirm new password** field do not exactly match the contents of the **New password** field. If this happens, correct the passcode confirmation.

Old password	<input type="text" value="MKLIZAMBVLFLBN"/>	
New password	<input type="text" value="MKLIZAMBVLFLBN"/>	valid
Confirm new password	<input type="text" value="MKLIZAMBVLXYZA"/>	unmatch ←

- Click the **Update** button in the lower right corner. Wait while the passcode is updated.



- When passcode updating is complete, the digitizer will be automatically reset. When resetting is complete, Passcode Editor will automatically close.




Troubleshooting


Digitizer Diagnostics, Firmware Download and Passcode Editor have built-in problem detection and troubleshooting help. After launching, each program first determines whether it can communicate with the digitizer. If communication is not possible, the program provides a checklist to help you identify and correct problems. As you complete each step, click its checkbox.

After each step, the program again attempts to communicate with the digitizer. If communication is established, you will be able to start normal program operations. If communication is not established, the next step in the troubleshooting process will be displayed.


1. Make sure the digitizer is connected to the computer with the appropriate USB or SCSI cable.



2. Turn off the digitizer, then turn it on.




3. Wait three seconds.





4. Scan for hardware changes:


a. Right-click **My Computer**, then click **Manage** in the pop-up menu.

- b. In the **Computer Management** window's left tree, click **Device Manager**.
- c. In the right tree, click the computer's name (at the top of the tree).
- d. In the toolbar, click **Scan for hardware changes** button.



5. In the **Computer Management** window, look for **Vidar Film Digitizer**, then click the corresponding checkbox in Digitizer Diagnostics.
6. To refresh this screen click on the icon the **Red Arrow** shown below.
7. You will find the digitizer listed either under Imaging Devices or Jungo depending on your operating system.





8. If you reach this point and do not find the Digitizer, no further troubleshooting suggestions are available. You can:
- Cancel the operation and Restart both the computer and digitizer (recommended).
 - Click the **Contact form** button in the lower right corner to send a report to Vidar Technical Support. Continue at step 7 in the “Digitizer Diagnostics” chapter of this manual. See Screen Shots below

? _ X

Digitizer Diagnostics 2.2

Model: _____
 Serial number: _____

System check

Requirements check

Run digitizer tasks

Done

Contact Information

Please fill out the form below and click the Next button below. Fields in bold are mandatory.

Digitizer Serial No.

Contact Name

Company Name

Address

Phone (Include country and area codes)

Email Address

Comments

<< BACK
NEXT >>



Digitizer SN: 340016

Film Digitizer Report - 2/22/2012 2:40:30 PM

NOTE: Please forward this report to VIDAR's Technical Support Group using one of the following options.

1. Email the report to medtech@vidar.com
2. Fax the report to the attention of Medtech at 703-471-7665.

Attn: Medical Technical Support
 365 Herndon Parkway
 Herndon, VA 20170
 +1.703.471.7070 (Phone)
 +1.703.471.7665 (Fax)
 MedTech@vidar.com



Digitizer Diagnostics 2.2

Model:
Serial number:

System check

Checking the digitizer's connection, please wait . . .

Unable to find a valid digitizer.

Exhausted digitizer connection troubleshooting options.

Troubleshooting Steps

- Check the digitizer USB or SCSI cable.
- Power the digitizer off, wait 3 seconds and then on.
- Wait for one minute.
- Try changing the hardware changes from the Device Manager.

Do you see a VIDAR Film Digitizer?

No Yes Try options 1 or 2 below.

1. Cycle power on the computer and digitizer and try again (recommended).
2. Click contact form button above to fill out contact form.

Send the report to MedTech for assistance.

CONTACT FORM

Film Digitizer Information

Hardware driver config	
Daver version	10.20
Toolkit version	
Model name	
Password	
Limited use device	
Maximum bit depth	
Maximum resolution	
Hardware version	
Kernel	
Firmware version	
DSP version	
MCB version	
Lamp type	
Lamp module	
Ballast time	
Feeder sn	
Feeder version	
Camera board	

Film Digitizer USB Information

Device description	Not found
Device bus speed	N/A

Appendix

This appendix reproduces a typical report generated by Digitizer Diagnostics.

Reports are initially created in the *My Documents\Film Digitizer\Report* folder. You can move or copy them to any other location. Note that if an email application is detected on the host pc the program will attempt to generate an email automatically for the report to be emailed to VIDAR. You do not have to use this feature and can cancel it if you desire.

Report names have the following syntax:

serial number – the number of times the report was run.pdf

For example, a report created for digitizer serial number 350282 would have this name:

350282 - 1.pdf

If additional reports are run the – number will change:

For example

350282 – 2.pdf

350282 – 3.pdf etc.