


# Scanner vs. Digitizer Images



VIDAR Dental Film Digitizer




Photo on Lightbox




In the scanner image the coping and crown are indiscernible; anterior mandibular teeth roots are washed out. Lose detail in lighter areas.



Photo on Lightbox



**VIDAR Dental**



**Scanner Image**

Lose detail in darker areas.

Scanner created streaks in image.



**Photo on Lightbox**

Diagnostic quality up to 3.6 out of 4.1 optical density.



**VIDAR Dental**




**Scanner Image**

Scanner can't capture higher densities.



Photo on Lightbox





Artifacts from dirty roller show with scanner, VIDAR digitizer nearly eliminates markings.